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
INSTRUCTIONS FOR USE OF ROTATING INSTRUMENTS AND SCREWDRIVERS

DRILLS: The dental drills "MODUS" are medical devices which have to be used as tools for drilling, flaring and/or tapping, after the evaluation and the planning of the implant surgery; The type of device, the diameters to be used, the lengths and the sequence of drilling (number of drills to be used) depend exclusively on the dentist, according to the surgical protocol that he want to follow; The maximum recommended speed is 500 rev/min.

- a) --- The circular drill (mucotome) is needed to create access to the bone surface as an alternative to dissection with a scalpel (flapless surgery); The maximum recommended speed is 300 rev/min;
- b)--- The initial step drill is used to cut the cortical region in a precise point for the subsequent drilling;
- c) --- The "standard" drills have reference marks, which identify the depth to reach; With these type of drills the dexterity and the experience of the surgeon are extremely important, especially for the stop at the selected and proper depth;
- d)--- The drills with "stop" have a mechanical stop (fixed or movable) that guarantees the depth of the drilling; With these drills we can avoid a deeper hole, which can be dangerous;
- e)--- The drills for guided implantology "Modus" must be guided by calibrated sleeves that ensure the direction and the depth of the insertion; These bushings eliminate the risk of drilling in wrong locations or even the risk of breaching the hole, with the consequent accidental invasion in cranial cavities and in nerve structures, that can cause long-lasting diseases, eg: temporary or permanent paresthesia of the semi-arch and of the innervation area, lesions to muscles and laceration of soft tissues; The presence on the drill body of a cylinder which has the same size for all the diameters and lengths, allows its use in guided implantology systems, where a simple and programmed sequence of drills is needed; The drills are designed in various diameters and various lengths; These drills can also be considered drills with "stop".
- The drilling diameters are designed in order to have a logical sequence that allows to get to the final drill, which must have the diameter (recommended by the manufacturer) of the implant to be inserted in the site.

COUNTER SINK: The counter sinks are used when there is the need to enlarge the initial part of the hole. They are made for adapting the shape to that of the neck of the implant to be inserted. The maximum recommended speed is 300 rev/min.

SCREW TAPS: In particularly dense bone (Type I) it is possible to use in advance a screw tap with the same profile of the implant to be inserted. The greater cutting power of the screw tap compared with that of the implant, allows you to prepare the site with less traumas. Normally it is recommended the manual tapping. If you prefer to use the dental micromotor, the maximum recommended speed is 30 rev/min.

TECHNICAL DATAS: All drills of the MODUS line are made of steel AISI 630 17-4 PH. They are subjected to an hardening heat treatment up to the H900 condition; They have a DLC coating.
Maximum number of device usage: 15

WARNINGS AND GENERAL PRECAUTIONS

- The non-observance of the surgical protocol that establishes diameters, lengths and sequence can cause serious damages to the patient, especially if you use drills of greater length;
- Check that the drills to use are in good condition, previously cleaned and sterilized;
- Ensure that drills are in good condition and no older than 15 uses;
- Before their use, check that the micromotor hand-piece fixes perfectly the drill and verify that the direction of rotation is the correct one;
- Ensure that the irrigation is adequate. It is required a copious irrigation with a sterile solution during drilling, in order to avoid damages to the bone tissue, with consequent bone necrosis;
- Do not exceed the maximum number of rev/min indicated on the label;
- The application of lever forces during drilling could cause breakage of the drilling or of the handpiece;
- When drilling always exercise an alternate pressure, using the technique of intermittent drilling;
- Always make sure that the laser writing indicating diameter and length is visible;
- The contingent eccentricity of the drill or its non-rectilinear shape may cause the over-sizing of the drilled hole;
- Always wear gloves when handling contaminated instruments and Wear eyes protection to screen yourself against particles.

SCREWDRIVERS AND INSTRUMENTS

The instruments manufactured by Errecieffe are reusable medical devices intended for transient use in the oral cavity. The functions of the screwdrivers is to tighten and unscrew all the connecting screws (cover screws, healing screws, screws for abutments, prosthetic screws, transfer screws, etc.). They are divided into:

- CONTRA-ANGLED SCREWDRIVERS FOR IMPLANTS AND PROSTHETIC: Suitable for the screwing with the handpiece;
- RATCHET SCREWDRIVERS FOR IMPLANTS AND PROSTHETIC: Indicated to continue screwing with the ratchet. Ratchet screwdrivers have a red polymer o'ring in the hexagon of the connection that ensures the friction between the tools;
- MANUAL SCREWDRIVERS FOR PROSTHETIC: Suitable to tightening the prosthetic screws manually;
- INSTRUMENTS (probe, ratchet, keys, drill extension shaft, parallel pins, tweezers).

WARNINGS: Excessive torques may strip the wells of the connecting screws and pare off the corners of the screwdrivers, causing even serious operative or prosthetic complications. The recommended torques for the various components are summed up as follow:

- Cover screws, healing screws, screws for MUA angled abutment: 15 Ncm
- Prosthetic screws: 25-32 Ncm
- Dental implants: 35-50 Ncm


Lever movements should be avoided as they increase the risk of breakage. Before tightening, make sure the screw head on the driver tip is correctly inserted into the screws to be tightened. Incorrect insertion is likely to pare off the connection of the screwdriver or the screw to be tightened. Replace drivers regularly.

These medical devices are delivered NON-STERILE. Before use, sterilize in steam autoclaves for 20 minutes at 134 ° C.

All rotating devices before being used must be cleaned, disinfected and sterilized. These processes are described in detail in our "Operating Procedure for cleaning, disinfection and sterilization of dental drills and endodontic instruments"

Storage and Handling: Devices should be adequately stored, protected from any dangerous situations.

Warning: the Italian legislation restricts the sale of this product only to dentists.

		<p>NON STERILE www.errechieffe.com E-mail: info@errechieffe.com</p>
---	---	--

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