-IT-

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UNIMAX[®]

by Biomax

INSTRUCTIONS FOR USE ZYGOMA AND PTERYGOID IMPLANT

MANUFACTURER IDENTIFICATION

The Manufacturer and Person responsible for placing on the market in the EU the Medical Devices called Dental Implants covered by this Instructions for Use is:

Errecieffe S.r.l.

Via V. Emanuele II, 68 - 24036 PONTE SAN PIETRO (BG) ITALY

Tel. +39 340 5181203 www.errecieffe.com info@errecieffe.com

DEVICE IDENTIFICATION

Zygomatic and Pterygoid dental implants of the ERRECIEFFE line with related surgical screws (where provided) and related implantable prosthetic components.

Intended use:

Dental implants (implant fixtures) are implantable devices designed to rehabilitate patients suffering from partial or total edentulism for the restoration of masticatory, phonetic or aesthetic function. Zygomatic implants are dental implants intended for the zygomatic bone. RCF zygomatic implants are suitable for immediate loading when good primary stability is achieved and with appropriate occlusal loading.

The pterygoid implant is intended for surgical placement in the upper jaw bone to provide support for prosthetic devices to restore chewing function. Pterygoid implants are only indicated for multiple restorations in splinted applications using at least two implants.

MATERIALS

The devices are made in:

Pterygoid: Titanium grade 4

Zygoma: Titanium grade 5

Instructions for use / Indications

Zygomatic titanium implants have been designed for all situations in which posterior mandibular atrophy complicates or prevents the use of conventional implants, representing a valid alternative to bone augmentation, sinus lift and bone transplantation in patients with atrophy in the posterior area of the maxilla. Zygomatic implants can solve the problem when bone is scarce, allowing the reconstruction of an entire arch with a minimum number of implants, and with a time saving of months.

The pterygoid implant is an endosseous fixture 18-28 mm long positioned in the region of the upper maxillary tuber; The longer length of the implant, compared to a standard fixture, allows you to reach the pterygoid fossa and therefore the sphenoid pillar on which the anchoring of the apical part of the implant takes place. The pterygoid implant is generally used in conjunction with the placement of other fixtures or teeth located in the mesial area with respect to the maxillary sinus.

DM Description

The ERRECIEFFE implant line is the alternative to implants already on the market, being "compatible" with others already on the market. The surgical procedures and warnings are the same as that normally used by the doctor.

Errecieffe has developed a surgical sequence for the insertion of the implants. This document, which specifies the type of drills to be used for site preparation, can be requested or downloaded from our website in the documents section.

The clinician must attend specific advanced training beyond the supervision of implantology experts on zygomatic implant placement.

RCF zygomatic implants are indicated to be surgically placed in the upper jaw, in the zygomatic bone, to support dental prostheses. They can be used to replace natural teeth in the event of loss due to both natural and traumatic causes.

RCF zygomatic implants can be tilted up to 60° with respect to the occlusal plane.

If used at angles between 30° and 60°, the inclined implants must be splinted; A minimum of 4 implants should be used when supporting a fixed denture in a fully edentulous arch.

The implant site begins with the exposure of the maxillary lateral wall, a mucoperiosteal flap, and the total thickness is lifted by making a ridge incision with bilateral distal vertical incision on the tuberosity areas. During surgical exposure of the lateral wall of the jaw it is absolutely necessary to pay attention to vital structures including nerves, veins and arteries. Injuries to these anatomical structures can lead to complications such as trauma to the eye, as well as extensive bleeding and dysfunction associated with nerves.

The placement of a pterygoid implant makes it possible to overcome the limitations dictated by the insertion of traditional implants in the upper molar regions characterized by limited residual volumes of the bone.

The posterior area of the tuberosity of the maxillary is closely connected at the bone level with the pterygoid process of the sphenoid and with the wall of the pyramidal process of the palatine, thus representing an appropriate anatomical structure, characterized by compact bone, suitable for supporting dental implants.

The pterygoid solution also makes it possible to avoid two-phase surgical procedures, i.e. a first phase of bone regeneration (sinus lift, bone grafts, split crest) followed, after 6-9 months, by the insertion of traditional dental implants; These procedures may not be accepted by the patient or may not give the desired results in terms of pre-implant regeneration.

The use of a pterygoid implant requires a high level of knowledge of the anatomy of the medical district of interest in order to design and execute the surgical act.

The RCF implant design has been thoroughly engineered with the help of specialized surgeons to facilitate insertion procedures, optimize primary stability, and maximize osseointegration.

CONTRAINDICATIONS

In evaluating the patient, in addition to considering the suitability for an implant-prosthetic rehabilitation, it is generally necessary to take into account the contraindications valid for dental surgery.

Absolute contraindications

The insertion of implants and implantable prostheses is contraindicated in patients with poor general health, poor or insufficient oral hygiene, impossibility or limited possibility of controlling their general conditions, or who have previously undergone organ transplants. Patients who are mentally ill, or who abuse alcohol, tobacco or drugs, who have serious neurological diseases, severe impairment of the immune system, chemotherapy in progress must also be discarded.

In the case of administration of bisphosphonates, numerous cases of peri-implant bone necrosis have been reported in the literature, mainly in the mandible. This problem particularly affects patients undergoing intravenous treatment.

Related contraindications

Patients with poor periodontal status must be previously treated and recovered or in the presence of infections and inflammations such as: periodontitis, gingivitis. In the event of a lack of bone substance or

poor quality of the recipient bone, such that the stability of the implant may be compromised, an appropriate guided tissue regeneration must be carried out beforehand. Metabolic or systemic diseases that compromise tissue regeneration. Intolerances or allergies to the material constituting the devices in question, dialysis, osteoporosis.

Pregnancy and postpartum period: you must wait for the end of pregnancy and the breastfeeding period. Follow-up tests should be carried out in the case of osteoporosis (blood values of acid and alkaline phosphatase).

WARNINGS

The risks of implant surgery include: perforation of the labial or lingual plate, bone fractures, implant fractures, fractures of the superstructures, cosmetic problems, inadvertent perforation of the nasal sinus, nerve injuries, compression of the natural dentition. The following pathophysiological problems may increase the risks: cardiovascular failure, chronic disorders, arrhythmia, chronic lung or respiratory diseases, gastrointestinal diseases, hepatitis, intestinal inflammation, chronic kidney failure, urinary system disorders, endocrine disorders, diabetes, thyroid disease, hematologic problems, anemia, leukemia, coagulation problems, osteoporosis or musculoskeletal arthritis, heart attack, neurological disorders, delays paralysis.

In case of damaged packaging, damaged or non-compliant device, do not use the device and notify the manufacturer or distributor of the problem.

Maximum importance must be given by users to the cleaning, disinfection and sterilization of the environment in which the intervention is carried out and of the instrumentation used to avoid the onset of infections or the occurrence of cross-infections

The devices in question are made of grade 4 titanium. Titanium is a non-ferromagnetic material that does not interact with the presence of magnetic fields. In the case of diagnostic investigations (e.g. MRIs) it is still advisable to notify the staff before performing the examination of the presence of devices implanted in Titanium. Do not use the devices in question on patients with intolerance or allergy to Titanium.

Users must carry out the surgery by correctly wearing all the protective devices required for their profession (e.g. mask, gowns, gloves, goggles) and pay the utmost attention to avoid small injuries from contact with devices with sharp or pointed portions and to avoid contamination with the patient's saliva or blood and to avoid contaminating the patient in turn.

Biomechanical considerations

Zygomatic and premaxillary fixtures, compared to standard implants, show a greater tendency to bend under transverse loads with respect to their axis. This is due to two factors: the significantly longer length and the fact that, in some cases, bone support is limited in the maxillary alveolar crest.

As a result, implants must be rigidly attached to standard implants fixed in the lower jaw. Based on clinical experience and biomechanical calculations, a total prosthesis of the jaw, supported by two zygomatic screws (one on each side), must be assisted by at least two regular stable implants in the anterior maxilla. The onset of bending moments, produced by forces transverse to the axis of the screws, is the most unfavorable condition and can potentially compromise the long-term stability of an implant-supported prosthesis. To reduce bending moments, the distribution of forces must be optimized through arch stabilization, reduced vestibular lever arms, reduced mesial/distal and anterior/posterior extensions, balanced occlusion, and reduced cusp tilt.

The dental implants in question are disposable devices supplied sterile. No reconditioning by the user is permitted. Once the sterile packaging has been opened, the implant must be used immediately, otherwise it must be thrown away.

Pay attention to the expiration date on the label of each package. Expired implants should not be used as sterility is not guaranteed.

No re-sterilization process is allowed by the user. Sterilization not carried out by the manufacturer does not guarantee the sterility of the device.

SIDE SYMPTOMS

Among the manifestations that accompany surgery, temporary complications such as pain, swelling, pronunciation problems, gingivitis, temporary local swelling, edema, hematomas, temporary limitations of sensitivity, temporary limitations of masticatory functions, post-operative micro haemorrhages in the following 12/24 hours may occur. After dental implant operations, the following may occur: loss of bone crest, permanent paresthesia, dysesthesia, local or systemic infections, exfoliation, hyperplasia, orofacial

and oronasal fistulas. As a precautionary measure, the patient should avoid activities that require physical exertion after surgery.

PREPARATORY TREATMENT INFORMATION

Proper planning of every aspect of oral surgery is essential for the long-term success of a prosthesis. The user of RCF zygomatic implant products has the duty to determine whether or not a product is suitable for the patient and the particular circumstances. Choosing inappropriate device sizes and implant malposition can lead to complications, side effects, and treatment failure. To avoid this, a detailed examination is necessary through clinical and radiographic analysis of the patient's general health condition and current medical status. Dental casts can be used to define the position and angle of implants. The dentist will develop the plan after collecting all the required data, because it provides important information and makes it possible to plan backwards, with the aim of improving the safety and outcomes of the surgical procedure.

The preparation phase for surgery with a ZYGOMATIC or PTERYGOID implant includes:

- Clinical picture and general medical and dental history
- Psychological and physiological assessment
- Examinations (e.g. radiological examination - CT scan)
- Evaluation of the general medical history, medical examinations, clinical examinations (complete blood count)
- Collection of information about the patient, with attention to factors that could interfere with the success of the treatment (e.g. cigarette smoking, poor oral hygiene, uncontrolled diabetes, active chemotherapy, bruxism)
- Professional dental cleaning
- Definition of an oral hygiene program, with possible periodontal interventions
- Provide information to patients about contraindications and unwanted side effects
- Prescription of antibiotic prophylaxis and possible pharmacological treatment, when necessary,
- Pre-prosthetic surgical planning and prosthetic planning in collaboration with the dental technician.
- Assessment of the risks associated with inadequate treatment of bone and mucosal tissues
- need for periodic post-checks

As with any implant patient case, an X-ray evaluation is essential. As far as the zygomatic implantation protocol is concerned, the presence of any pathology in the maxillary sinuses must be detected, bearing in mind that the thickness of the antral mucosa must not exceed 6 mm. It is also mandatory to check the patency of the ostiomeatal complex. The volume of the zygomatic bone must then be evaluated (the minimum recommended thickness is 6 mm).

SURGICAL INTERVENTION

Surgical techniques for implants are taught in universities to dental graduates. However, the following factors should be taken into account:

- tissues, both hard and soft, must be treated with extreme care, taking all the necessary precautions to obtain a good integration of the implant;
- the normal biological principles of osseointegration must be respected;
- thermal trauma must be avoided as they would necrotize and compromise the possibility of osseointegration. For this purpose, adequate drilling speeds must be used, drills with a cutting edge in excellent condition, drilling must be carried out intermittently, cooling the site with the necessary irrigation and the hole must be enlarged using drills with progressively larger specific diameters.
- It is essential to respect the healing times recommended in implant surgery and to periodically check, even with radiographic checks, the progressive state of osseointegration.

Everything possible must be done to minimize damage to the host tissue. Particular attention must be paid to thermal and surgical trauma and the elimination of contaminants and sources of infection. The surgical procedure requires a high degree of precision and care. Any deviation from the principle of as little trauma as possible during implant installation increases the risk of failure to achieve osseointegration. The use of sharp drills, sufficient irrigation, a drilling motion in and out, short cutting cycles, waiting for the bone to cool down, and the use of pilot drills of increasing size are essential.

Special care should be taken to avoid over- or under-preparation of the osteotomy. The implants must be inserted in such a way that they are stable and free of mobility. Excessive insertion torque (greater than

60 Ncm) can cause damage to the implant or instruments and fracture or necrosis of the bone site. All instruments used in surgery must be kept in good condition and care must be taken that the instruments do not damage the implants or other components. Precautions should be taken to avoid ingesting or aspirating components used in implantology. After the implant is installed, the surgeon's assessment of bone quality and initial stability will determine when it will be possible to load the implants. An appropriate follow-up protocol should be followed.

The surgical intervention for the zygomatic involves the following phases:

- Incision the ridge of the edentulous maxilla with distal vertical release incision.
- Reflect a full-thickness mucoperiosteal flap exposing the lateral maxillary wall.
- Expose the alveolar ridge, including the palatal side.
- Carefully dissect down to the level of the infraorbital foramen.
- Reflect laterally at the infraorbital nerve level and expose the body of the zygomatic bone.
- Caution: It is essential to identify and protect the infraorbital nerve.
- Place a retractor in the frontozygomatic notch to facilitate visualization of the intended apical point of the implant
- Make a window of about 10 mm x 5 mm on the lateral wall of the breast, near the infrazygomatic ridge.
- Lift the sinus mucosa from the area where the implant will pass by following the extra-sinus path, from the sinus floor to the roof, trying not to penetrate the membrane. However, penetration of the sinus membrane will not result in a negative outcome.
- Identify the implant trajectory considering that the prosthetic connection level of the zygomatic implant must emerge at the position of the second premolar near the edge of the ridge. The zygomatic implant must penetrate the base of the zygomatic bone as posteriorly as possible.
- To prepare the direction of drilling of the zygomatic bone, it is necessary to make an orientation channel in the jaw using diamond-coated drills. Maximum speed 1500 rpm.

The classic technique for zygomatic placement involves cutting a window of the maxillary sinus and placing the implant through the breast. Since then, the maxillary sinus slit technique and the exteriorized technique have been developed, with the implant placed through the sinus wall and outside the sinus wall, respectively. It has been suggested that the choice of technique should consider the ridge, ridge concavity, and sinus anatomy.

Anatomically, patients vary with regard to the concavity of the lateral wall of the breast and the zygoma in relation to the crest. The classification of zygomatic implants according to Aparicio presents the five different types according to the concavity of the jaw wall. In each type, the implant is placed on the crest inferiorly and on the zygoma superiorly. A part of the implant body may not be enclosed in the jaw wall, which means that, in some clinical situations, it is located outside the lateral wall of the breast.

Surgery for pterygoid involves the following phases:

Make a full-thickness ridge incision of the maxilla tuber and then vestibular and palatal release incisions. Proceed with the skeletonization of the posterior upper jaw and the identification of clinical and radiological anatomical findings.

Use the initial drill to create an invitation on the cortical bone useful for positioning the tip of the depth drill.

Use the depth drill to create the depth pilot hole, up to the pterygoid fossa.

Proceed to enlarge the hole using the correct drill.

DEVICE ACCESSORIES

Description of other devices used in conjunction with the device

Surgical screw (or screw cap)	It can be used in the post-operative phase and whose duration of use can normally be used for a continuous duration of more than 30 days. It is the Clinician's faculty to use it in the post-operative phase.
----------------------------------	--

Mounter	Used only as a transport device for inserting the fixture into the oral cavity, the duration of which can normally be used for a continuous duration of less than 60 minutes.
---------	---

Additional devices, which can be used with dental implants (covered in the relevant Technical Files), are:

Healing (or transmucosal) abutments	Small direct-screwing abutments, of different heights according to the different implant systems, intended to recondition the emergence profiles of the mucous membranes, before prosthetic loading. They are surgical invasive devices intended for the oral cavity, for long term.
Temporary abutments	Temporary abutments, normally composed of a titanium base with an upper cannula on which the dentist or dental technician re-bases an acrylic-type prosthesis. In some versions they are made in Peek, which can be modified by milling from the laboratory or by the doctor directly to the chair. Peek is not resin-based, so these abutments are normally used for the rehabilitation of single prostheses by cementing a crown.
Abutments for traditional screw-retained prostheses	Straight and angled abutments for cemented prostheses. Millable abutments, preformed for cemented prostheses. They are customized by milling by the dental laboratory or by CAD-CAM techniques by milling centers. Castable pillars with preformed alloy base. They are used to obtain individual abutments for cemented prostheses or for overdenture bars or structures for screw-retained bridges such as Toronto Bridge by overcasting in the dental laboratory. Titanium bases equipped with connection to the implant, which in the upper part have a standard coupling cone for CAD-CAM techniques. Abutments that are preferably used for screwing multiple prostheses (Toronto-type prostheses) of the traditional type (straight with direct screwing) and for the MUA technique (for the correction of large disparellelisms, complete with the relative screws for fixing them on the implants).
Removable Overdenture Attachment Components	Spherical attachments that function as buttons for the stabilization of a total prosthesis. Spherical attachments require that a suitable matrix is positioned inside the prosthesis in correspondence with the attachments capable of hooking onto the spherical heads of the attachments themselves. These males are made up of caps in polyamide, or in gold alloy, or in titanium or by O-Ring type attachments.
Transfer	They have the function of transferring the exact position of the implant connection from the mouth to the dental model, in terms of height, inclination and indexing. There are various types of transfers: tear-off, pick-up technique transfer and (limited to some implantological systems) pull-up type. Special transfers (CAD-CAM transfers) are also available, intended to be screwed onto the plaster models and not in the mouth, which allow the use of three-dimensional software to transfer the position of the platforms from the plaster models to virtual CAD models for the subsequent techniques for making individual prostheses using the CAM process.
Castable abutments	They are used to obtain, by casting in the dental laboratory, individual pillars for cemented prostheses or for the fusion of bars for overdentures or structures for screw-retained bridges such as Toronto Bridge.
Tightening screws for abutments and Superstructures	These are the screws needed to screw pillars and superstructures.

USERS

The use and handling of the medical device is reserved for medical and dental personnel with the necessary professional qualification and training.

THE CHOICE OF THE DEVICE MUST BE CONSEQUENT TO AN ACCURATE PATIENT HISTORY

Specific responsibilities

Errecieffe products for oral implantology meet the GENERAL SAFETY AND PERFORMANCE REQUIREMENTS established by Regulation (EU) 2017/745-MDR. Any use of the products other than the specific one is to be considered as "improper use", relieving the manufacturer of any liability.

DISPOSAL PROCEDURES

The components, if removed from the oral cavity due to biological or mechanical failure, must be assimilated to biological waste, according to the regulations in force at local level.

INCIDENT REPORTING

In accordance with SECTION 2 VIGILANCE, Article 87 Reporting of serious incidents and safety corrective actions of MDR 745/2017/EU, DL 15 November 2005 and MEDDEV guideline 2.12-1, it is recalled that public or private healthcare professionals are required to report accidents or near misses with Medical Devices to the manufacturer (Errecieffe S.r.l.) or to the National Competent Authorities as a matter of urgency.

An accident is defined as any dysfunction or deterioration in the characteristics and/or performance of a device, as well as any deficiency in labelling or instructions for use which, directly or indirectly, may cause or have caused death or a serious deterioration in the state of health of the patient or a user or other persons.

IMPLANT PASSPORT

Errecieffe S.r.l. provides, together with the Medical Device, the card for the dental implant with the related information that is provided to the patients.

SUMMARY DOCUMENT ON SAFETY AND CLINICAL PERFORMANCE

See SUMMARY document at EUDAMED website

FURTHER INFORMATION

This device does not contain software.

This device must be sterilized before use (the device is supplied in NON-STERILE condition).

This device is intended for re-use.

This device is not intended for use by lay users.

The device is not listed as a product that does not have a medical use.

LIABILITY FOR DEFECTIVE PRODUCT AND WARRANTY TERMS

Optimal patient care and attention to his needs are necessary conditions for implantological success and it is therefore necessary to carefully select the patient, inform him of the inherent risks and duties associated with the treatment and encourage him to cooperate with the dentist for the success of the treatment itself. It is therefore necessary that the patient maintains good hygiene, confirmed during check-ups and follow-up appointments; it must always be ensured and documented as, moreover, the indications and prescriptions of the doctor before and after surgery must be observed and documented.

The instructions provided by Errecieffe S.r.l. are available at the time of treatment and accepted by dental practice; they must be observed and applied at all stages of care: from the patient's medical history to post-operative check-ups.

The warranty covers only ascertained production defects, upon shipment of the part identified by article number and batch, within the warranty period. The warranty clauses are available on the www.errecieffe.com website.

This device must not be sterilized before use (the device is supplied in a STERILE condition).

This device is not intended to be reused (the device is DISPOSABLE). This device is not intended for use by laypeople.

The device is not included in the list of product groups that do not have a medical use.

LEGEND OF SYMBOLS USED ON PACKAGING

SYMBOL	DESCRIPTION	SYMBOL	DESCRIPTION
	Manufacturer's Name		CE marking with identification of the Notified Body
	Sterile by irradiation		Non-sterilizable product
	Single-use product, do not reuse		Do not use if the packaging is damaged
	Consult the instructions for use		Caution! See instruction for use
	Device ID - Code		Batch number
	Unique Device Identification - Alphanumeric code that uniquely identifies a medical device		The product is a Medical Device
	Expiration Date		Distributor

DATE AND VALIDITY OF THESE OPERATING INSTRUCTIONS

These operating instructions are valid and effective from February 2025.

BIOMAX SPA

VIA ZAMENHOF 615 - 36100 VICENZA
T.0444913410 - info@biomax.it